

BIOMEDICAL RESEARCH SERVICE CENTER

UNIVERSITY at BUFFALO, STATE UNIVERSITY of NEW YORK

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 3435 Main Street, Buffalo, NY 14214, USA
Tel/Fax: (716) 829-3106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

Aspartate Aminotransferase (AST/GOT) Assay Kit (Cat #: E-116)

COMPONENTS: AST Substrate- 5 ml, store at -70°C
TA Assay Solution- 10 ml (for 200 wells), store at -70°C
10x Cell Lysis Solution- 25 ml, store at 4°C (**contains 1% TX-100; swirl bottle briefly prior to dilution**)

PRODUCT DESCRIPTION: The assay is based on sequential AST-mediated transamination reaction and glutamate dehydrogenase reaction, which couples the reduction of INT to formazan ($\epsilon = 18 \text{ mM}^{-1}\text{cm}^{-1}$) for detection of AST activity in cell/tissue extracts and serum/plasma. Kit components are stable for several years if stored and handled properly.

Preparation of cell/tissue extracts:

1. Prepare 1x Cell Lysis Solution by diluting 10x Cell Lysis Solution with ice-cold dH₂O. Bring up at least $\sim 10^5$ washed cells in 50 – 100 μl ice-cold 1x Cell Lysis Solution by pipetting up and down gently. Leave lysate on ice for 5 min with agitation. If lysate is overly turbid, add more 1x Cell Lysis Solution and repeat pipetting. Tissue is homogenized in ice-cold 1x Cell Lysis Solution ($\sim 10 \text{ mg}$ tissue in 0.5 ml).
2. Centrifuge lysate in a cold microfuge at $\sim 14,000 \text{ rpm}$ for 5 min. Supernatant is harvested and stored at -80°C.
3. Use the BCA protein assay method to determine lysate protein concentration. Dilute sample protein concentration to 0.1 mg/ml using ice-cold 1x Cell Lysis Solution.

Dilution of serum/plasma:

Each serum/plasma sample should be diluted 10-fold with ice-cold 1x Cell Lysis Solution prior to assay (e.g. 90 μl 1x Cell Lysis Solution + 10 μl sample). Mix well. Use 5 μl diluted sample for AST assay. Do not use undiluted serum/plasma for the assay.

Reagent thawing:

Keep thawed AST Substrate and TA Assay Solution on ice (do not over thaw). Vortex substrate solution to re-dissolve contents if necessary. Minimize the time that the reagents are thawed. Freeze solutions immediately after use.

Enzyme activity assay:

1. Each sample to be assayed has a control well (containing 5 μl sample and 45 μl dH₂O) and a reaction well (containing 5 μl sample and 45 μl AST Substrate). Pipette 5 μl of each sample to a 96-well plate in duplicate: one set for control wells and another set for reaction wells.
2. After all samples have been pipetted to the plate in duplicate, add 45 μl dH₂O to each control well and 45 μl AST Substrate to each reaction well. Gently agitate plate for 10 sec to mix. Cover plate and incubate in a 37°C incubator for 60 min. Do not use CO₂ incubator. This incubation step allows transamination reaction to take place.
3. Remove plate from the incubator. Add 50 μl TA Assay Solution to each control and reaction well. Gently agitate plate 10 sec to mix. Cover plate and incubate in a 37°C incubator for another 60 min. This step measures the amount of glutamate generated by AST.
4. Remove plate from the incubator. Terminate reaction by adding 50 μl 3% Acetic acid (not included in the kit) to each control and reaction well followed by brief gentle agitation. Measure O.D._{492 nm} using a plate reader. Subtract control well reading from reaction well reading for each sample. Use the subtracted sample reading ($\Delta\text{O.D.}$) for enzyme activity calculation as shown below.
5. AST enzyme activity in IU/L unit = $\mu\text{mol}/(\text{L}\cdot\text{min}) = (\Delta\text{O.D.} \times 1000 \times 150 \mu\text{l}) \div (60 \text{ min} \times 0.7 \text{ cm} \times 18 \times 5 \mu\text{l})$
= $\Delta\text{O.D.} \times 39.68$. For tissue/cell extracts, enzyme activity can be presented as units/ μg proteins.

Additional information:

- A solution of 3% Acetic acid needs to be prepared for reaction termination.
- The assay solution contains DMSO and iodinitrotetrazolium violet. Please refer to the product page of our website or contact us for MSDS information.