

BIOMEDICAL RESEARCH SERVICE CENTER

UNIVERSITY at BUFFALO, STATE UNIVERSITY of NEW YORK

Department of Biochemistry, Attn: Dr. Lee, University at Buffalo, 3435 Main Street, Buffalo, NY 14214, USA
Tel/Fax: (716) 829-3106 Email: chunglee@buffalo.edu Web: www.bmrservice.com

NAD⁺/NADH Assay kit (Cat #: A-118)

COMPONENTS: NAD⁺/NADH Assay Solution- 5 ml (for 100 wells), store at -70°C (**avoid light exposure**)
NADH Assay Solution- 5 ml (for 100 wells), store at -70°C (**avoid light exposure**)
NADH- 20 mg, store at -70°C (**avoid light exposure**)
NADH Buffer- 1.5 ml, store at 4°C
PEG Solution: 5 ml; store at 4°C (**viscous; pipette solution with a cut tip**)

PRODUCT DESCRIPTION: Nicotinamide adenine dinucleotide exists in two forms in cells: NAD⁺ and NADH, which are involved in numerous redox reactions. The NAD⁺/NADH Assay Kit is based on LDH-mediated reduction of NAD⁺ and the tetrazolium salt INT in a diaphorase-coupled reaction, which generates a formazan product that can be measured by absorbance at 492 nm (detection limit 1 – 5 µM). Reagents are stable for at least 1 year if stored and handled properly.

PROTOCOLS

Preparation of tissue/cell samples: Cell and tissue samples are deproteinized by PEG precipitation. Please follow the extraction protocol at <http://www.bmrservice.com/SupplementPEG.html>. Alternatively, the samples can be deproteinized by TCA precipitation (<http://www.bmrservice.com/SupplementTCA.html>), which is recommended for analysis of nucleotides. Some TCA-treated samples may be diluted with dH₂O prior to assay to achieve optimal result.

Reagent thawing:


Keep thawed Assay Solutions on ice and shielded from light during assay. Do not over thaw solutions. It is important to minimize the time that the reagents are thawed. Freeze reagents immediately after use.

NADH standards:

Make fresh a 20 mM NADH standard: carefully weigh 1 mg of NADH and vortex to dissolve NADH in 70 µl of ice-cold NADH Buffer in a tube shielded from light. Dilute the standard 200-fold with NADH Buffer to generate 100 µM NADH, e.g., 199 µl NADH Buffer + 1 µl 20 mM NADH. Perform additional 1:1 dilution in NADH Buffer to generate 50 µM, 25 µM, and 12.5 µM NADH. Keep NADH solutions on ice shielded from light during assay.

ASSAY

1. Use a plain (uncoated) 96-well plate for the assay. Add 20 µl of dH₂O to a well. Add 20 µl of each NADH standard to a set of wells. Add each deproteinized sample to two sets of wells (**in duplicate**).
2. Gently agitate Assay Solutions prior to first pipetting. Add 50 µl NADH Assay Solution to each NADH standard well and to one set of sample wells. Add 50 µl NAD⁺/NADH Assay Solution to the other set of sample wells and the dH₂O well. Mix contents by gentle but thorough agitation for 10 sec.
3. Cover and incubate plate in a 37°C humidified incubator for 60 – 90 min. Do not use CO₂ incubator.
4. Measure absorbance at 492 nm using a plate reader. Eliminate air bubbles present in wells prior to measurement. Note: O.D. will be inflated if turbidity appears in wells. Samples can be diluted with dH₂O prior to assay to eliminate turbidity. Some TCA-treated samples may be diluted with dH₂O prior to assay to achieve optimal result.
5. Subtract the dH₂O well reading from all readings of the NAD⁺/NADH wells prior to calculation (see below).
6. Plot NADH standards vs. O.D._{492 nm}. Generate a trendline equation on chart. Use the equation to convert sample O.D. from the two sets of measurements to sample [NADH] and [NADH + NAD⁺]. x = concentration in µM; y = O.D._{492 nm}. A new plot must be generated for each assay.
7. Sample NADH concentration is derived from the NADH Assay Solution wells. Sample (NADH + NAD⁺) concentration is derived from the NAD⁺/NADH Assay Solution wells. Subtraction of NADH concentration from (NADH + NAD⁺) concentration yields sample NAD⁺ concentration in µM.



Additional Information:

- Assay Solutions contain the organic solvent DMSO and idonitrotetrazolium violet. Please contact us or visit the product webpage for MSDS information.